

Anti-Hoxc4 Polyclonal Antibody For Wb, Elisa - P09017

Item Number: CM0015531

Introduction

Rabbit polyclonal antibody against human HOXC4, a developmental transcription factor for anterior-posterior patterning. Validated in WB and ELISA, with cross-reactivity to mouse and rat. Ideal for developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	HOXC4 Rabbit pAb
Remarks/Alias	HOX3, cp19, HOX3E, HOXC4
Species	Human
GeneID (Human)	3221
GeneID	3221
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-130 of human HOXC4 (NP_055435.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P09017
Protein Weight	30kDa
Shipping	Ice bag