

Anti-Hoxc6 Rabbit Monoclonal Antibody For Wb, Elisa - P09630

Item Number: CM0007601

Introduction

Rabbit monoclonal antibody targeting human HOXC6, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying developmental transcription regulation and anterior-posterior axis patterning.

[Learn More](#)

Parameter	Value
Product Name	HoxC6 Rabbit mAb
Remarks/Alias	CP25; HOX3; HOX3C; HHO.C8; HoxC6
Species	Human
Gene ID (Human)	3223
Gene ID	3223
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human HoxC6 (NP_004494.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09630
Protein Weight	27kDa
Shipping	Ice bag