

Anti-Hoxd10 Polyclonal Antibody For Wb, Elisa - P28358

Item Number: CM0030580

Introduction

Rabbit polyclonal antibody against human HOXD10, validated for WB and ELISA. Detects endogenous levels of HOXD10. Suitable for human, mouse, and rat samples. Ideal for developmental biology and urogenital tract research.

[Learn More](#)

Parameter	Value
Product Name	HOXD10 Rabbit pAb
Remarks/Alias	HOX4; HOX4D; HOX4E; Hox-4.4; HOXD10
Species	Human
Gene ID (Human)	3236
Gene ID	3236
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human HOXD10 (NP_002139.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28358
Protein Weight	38kDa
Shipping	Ice bag