

Anti-Hoxd10 Rabbit Monoclonal Antibody For Wb, Elisa - P28358

Item Number: CM0006988

Introduction

Anti-HOXD10 rabbit monoclonal antibody for human HOXD10. Validated for WB and ELISA, cross-reacts with mouse. Ideal for developmental biology and urogenital tract research.

[Learn More](#)

Parameter	Value
Product Name	HOXD10 Rabbit mAb
Remarks/Alias	HOX4; HOX4D; HOX4E; Hox-4.4; HOXD10
Species	Human
GeneID (Human)	3236
GeneID	3236
Immunogen	Recombinant protein corresponding to a sequence within amino acids 1-150 of human HOXD10 (NP_002139.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P28358
Protein Weight	38kDa
Shipping	Ice bag