

Anti-Hoxd8 Polyclonal Antibody For Wb, Elisa - P13378

Item Number: CM0022688

Introduction

Rabbit polyclonal antibody targeting HOXD8, a homeobox transcription factor crucial for anterior-posterior patterning. Validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	HOXD8 Rabbit pAb
Remarks/Alias	HOX4; HOX4E; HOX5.4; HOXD8
Species	Human
GeneID (human)	3234
GeneID	3234
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200 to the C-terminus of human HOXD8 (NP_001186675.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13378
Protein Weight	32kDa
Shipping	Ice bag