

Anti-Hoxd9 Polyclonal Antibody For Wb, Elisa - P28356

Item Number: CM0019871

Introduction

Rabbit polyclonal antibody targeting human HOXD9, a homeobox transcription factor critical for developmental anterior-posterior axis specification. Validated for Western blot (WB) and ELISA applications. Detects a 36kDa band, cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	HOXD9 Rabbit pAb
Remarks/Alias	HOX4, HOX4C, Hox-4.3, Hox-5.2, HOXD9
Species	Human
GeneID (Human)	3235
GeneID	3235
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human HOXD9 (NP_055028.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P28356
Protein Weight	36kDa
Shipping	Ice bag