


Anti-Hp1 Alpha/Cbx5 Polyclonal Antibody For Wb, If/Icc, Elisa - P45973

Item Number: CM0015781



Introduction

High-quality anti-HP1 alpha/CBX5 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Recognizes human and mouse HP1 alpha, a heterochromatin protein 1 homolog involved in epigenetic repression and kinetochore formation. Ideal for chromatin research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HP1 alpha/CBX5 Rabbit pAb
Remarks/Alias	HP1; HP1A; HEL25; X5
Species	Human
Gene ID (Human)	23468
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human HP1 alpha/CBX5 (NP_001120794.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P45973
Protein Weight	22kDa
Shipping	Ice bag