


Anti-Hp1 Alpha/Cbx5 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P45973

Item Number: CM0010417



Introduction

Human HP1 alpha/CBX5 rabbit monoclonal antibody for WB, IHC-P, IF/ICC, IP, ELISA. Detects human, mouse, rat targets. Suitable for heterochromatin and epigenetic research. P45973.

[Learn More](#)

Parameter	Specification
Product Name	HP1 alpha/CBX5 Rabbit mAb
Remarks/Alias	HP1; HP1A; HEL25; HP1 alpha/CBX5
Species	Human
GeneID (Human)	23468
GeneID	23468
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 92-191 of human HP1 alpha/CBX5 (P45973).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P45973
Protein Weight	22kDa
Shipping	Ice bag