


Anti-Hp1 Alpha/Cbx5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa, Chip - P45973

Item Number: CM0018072



Introduction

Rabbit polyclonal antibody targeting HP1 alpha/CBX5, applicable for WB, IHC-P, IF/ICC, IP, ELISA, and ChIP. Recognizes human, mouse, and rat HP1 alpha, a heterochromatin protein involved in epigenetic repression. Ideal for chromatin and centromere research.

[Learn More](#)

Parameter	Value
Product Name	HP1 alpha/CBX5 Rabbit pAb
Remarks/Alias	HP1; HP1A; HEL25; HP1 alpha/CBX5
Species	Human
GeneID (Human)	23468
GeneID	23468
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-191 of human HP1 alpha/CBX5 (NP_001120794.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P45973
Protein Weight	22kDa
Shipping	Ice bag