

Anti-Hp1 Gamma/Cbx3 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13185

Item Number: CM0009698



Introduction

Rabbit monoclonal antibody against human HP1 gamma/CBX3, validated for WB, IHC-P, ELISA in human, mouse, rat. Suitable for chromatin and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	HP1 gamma/CBX3 Rabbit mAb
Remarks/Alias	HECH; HP1-GAMMA; HP1Hs-gamma; HP1 gamma/CBX3
Species	Human
GeneID (human)	11335
GeneID	11335
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human HP1 gamma/CBX3 (NP_057671.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13185
Protein Weight	21kDa
Shipping	Ice bag