

Anti-Hp1 Gamma/Cbx3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13185

Item Number: CM0024008



Introduction

High-quality rabbit polyclonal antibody targeting human HP1 gamma/CBX3 (HECH, HP1-GAMMA) for reliable detection in WB, IHC-P, and ELISA across human, mouse, and rat samples. Validated for epigenetics and heterochromatin research.

[Learn More](#)

Parameter	Value
Product Name	HP1 gamma/CBX3 Rabbit pAb
Remarks/Alias	HECH; HP1-GAMMA; HP1Hs-gamma; HP1 gamma/CBX3
Species	Human
GeneID	11335
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-183 of human HP1 gamma/CBX3 (NP_009207.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13185
Protein Weight	21kDa
Shipping	Ice bag