

Anti-Hp1 Gamma/Cbx3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q13185

Item Number: CM0023189



Introduction

Rabbit polyclonal antibody targeting HP1 gamma/CBX3, validated for WB, IHC-P, IF/ICC, IP, ELISA. Recognizes human, mouse, rat. Molecular weight: 21kDa. Ideal for heterochromatin and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] HP1 gamma/CBX3 Rabbit pAb
Remarks/Alias	HECH; HP1-GAMMA; HP1Hs-gamma; X3
Species	Human
GeneID (Human)	11335
GeneID	11335
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-183 of human HP1 gamma/CBX3 (NP_009207.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13185
Protein Weight	21kDa
Shipping	Ice bag