

Anti-Hp1Bp3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q5Ssj5

Item Number: CM0014687



Introduction

HP1BP3 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets heterochromatin protein 1-binding protein 3, a regulator of cell proliferation and chromatin condensation.

[Learn More](#)

Parameter	Value
Product Name	HP1BP3 Rabbit pAb
Remarks/Alias	HP1BP74; HP1-BP74; HP1BP3
Species	Human
GeneID	50809
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human HP1BP3 (NP_057371.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5SSJ5
Protein Weight	61kDa
Shipping	Ice bag