


## Anti-Hpcal1 Polyclonal Antibody For Wb, Elisa - P37235

Item Number: CM0018472



### Introduction

HPCAL1 Rabbit pAb (P37235) is a polyclonal antibody validated for WB and ELISA, reactive with human, mouse, and rat. Ideal for calcium-binding protein and visual signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | HPCAL1 Rabbit pAb                                                                                                  |
| Remarks/Alias    | BDR1; HLP2; VILIP-3; HPCAL1                                                                                        |
| Species          | Human                                                                                                              |
| Gene ID (Human)  | 3241                                                                                                               |
| Gene ID          | 3241                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-193 of human HPCAL1 (NP_002140.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | P37235                                                                                                             |
| Protein Weight   | 22kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |