

Anti-HPD Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P32754

Item Number: CM0030725



Introduction

Rabbit polyclonal antibody against human HPD (4-hydroxyphenylpyruvate dioxygenase) for WB, IHC-P, ELISA. Confirmed cross-reactivity with mouse and rat. Useful in tyrosine metabolism and metabolic disorder research.

[Learn More](#)

Parameter	Value
Product Name	HPD Rabbit pAb
Remarks/Alias	PPD; 4HPPD; GLOD3; 4-HPPD; HPPDASE; HPD
Species	Human
Gene ID (Human)	3242
Gene ID	3242
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-393 of human HPD (NP_002141.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32754
Protein Weight	45kDa
Shipping	Ice bag