


Anti-Hprt1 Polyclonal Antibody For Wb, Elisa - P00492

Item Number: CM0018828



Introduction

Rabbit polyclonal antibody against human HPRT1, suitable for Western blotting and ELISA. Shows cross-reactivity with mouse and rat. Useful for studying purine salvage pathway, Lesch-Nyhan syndrome, and gout.

[Learn More](#)

Parameter	Value
Product Name	HPRT1 Rabbit pAb
Remarks/Alias	HPRT; HGPRT; HPRT1
Species	Human
GeneID (Human)	3251
GeneID	3251
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-218 of human HPRT1 (NP_000185.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00492
Protein Weight	25kDa
Shipping	Ice bag