


Anti-Hps4 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nqg7

Item Number: CM0019869



Introduction

Rabbit polyclonal antibody to HPS4, a component of the BLOC-3 complex. Validated for WB, IF/ICC, ELISA in human, mouse, rat. Ideal for studying melanosome biogenesis and Hermansky-Pudlak syndrome.

[Learn More](#)

Parameter	Value
Product Name	HPS4 Rabbit pAb
Remarks/Alias	LE; BLOC3S2; HPS4
Species	Human
GeneID (Human)	89781
GeneID	89781
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 387-708 of human HPS4 (NP_071364.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQG7
Protein Weight	77kDa
Shipping	Ice bag