

Anti-Hr Polyclonal Antibody For Wb And Elisa - 043593

Item Number: CM0013745

Introduction

Recombinant protein immunogen corresponding to amino acids 200-480 of human HR. Rabbit anti-HR polyclonal antibody validated for WB and ELISA; cross-reacts with human, mouse, rat. Suitable for studying histone demethylation and hair biology.

[Learn More](#)

Parameter	Value
Product Name	HR Rabbit pAb
Remarks/Alias	AU; MUHH; ALUNC; HYPT4; MUHH1; HSA277165; HR
Species	Human
GeneID (Human)	55806
GeneID	55806
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-480 of human HR (NP_005135.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43593
Protein Weight	127kDa
Shipping	Ice bag