


## Anti-Hras Monoclonal Antibody For Wb, Elisa - P01112

Item Number: CM0002508



### Introduction

Anti-HRAS rabbit monoclonal antibody validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying Ras-MAPK signaling, cell proliferation, and cancer biology. SWISS P01112.

[Learn More](#)

| Parameter        | Value                                                                            |
|------------------|----------------------------------------------------------------------------------|
| Product Name     | HRAS Rabbit mAb                                                                  |
| Remarks/Alias    | CTLO; HAMSV; HRAS1; RASH1; p21ras; C-H-RAS; H-RASIDX; C-BAS/HAS; C-HA-RAS1; HRAS |
| Species          | Human                                                                            |
| GeneID (Human)   | 3265                                                                             |
| GeneID           | 3265                                                                             |
| Immunogen        | Synthetic Peptide                                                                |
| Source           | Rabbit                                                                           |
| Category         | Monoclonal Antibodies                                                            |
| Application      | WB, ELISA                                                                        |
| Cross Reactivity | Mouse, Rat                                                                       |
| SWISS            | P01112                                                                           |
| Protein Weight   | 21kDa                                                                            |
| Shipping         | Ice bag                                                                          |