


Anti-Hras Polyclonal Antibody For Wb, If/Icc, Elisa - P01112

Item Number: CM0017052



Introduction

Rabbit polyclonal antibody targeting human HRAS (GTPase HRas, p21ras). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse. Suitable for studying Ras-MAPK signaling and oncogene research.

[Learn More](#)

Parameter	Value
Product Name	HRAS Rabbit pAb
Remarks/Alias	CTLO; HAMSV; HRAS1; RASH1; p21ras; C-H-RAS; H-RASIDX; C-BAS/HAS; C-HA-RAS1; HRAS
Species	Human
Gene ID	3265
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-189 of human HRAS (NP_005334.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01112
Protein Weight	21kDa
Shipping	Ice bag