


Anti-Hras Polyclonal Antibody For Wb, Ihc-P, Elisa - P01112

Item Number: CM0029232



Introduction

High-quality rabbit polyclonal antibody targeting human HRAS (p21ras). Validated for WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Ideal for Ras-MAPK signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	HRAS Rabbit pAb
Remarks/Alias	CTLO; HAMSIV; HRAS1; RASH1; p21ras; C-H-RAS; H-RASIDX; C-BAS/HAS; C-HA-RAS1; HRAS
Species	Human
GeneID (human)	3265
GeneID	3265
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human HRAS (NP_203524.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01112
Protein Weight	21kDa
Shipping	Ice bag