

Anti-Hrh1 Polyclonal Antibody For Wb, Elisa - P35367

Item Number: CM0013014



Introduction

Rabbit polyclonal antibody targeting human histamine H1 receptor (HRH1). Suitable for WB and ELISA applications. Detects a 56kDa protein. Ideal for GPCR signaling, allergy, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	HRH1 Rabbit pAb
Remarks/Alias	H1R; H1-R; HH1R; hisH1; HRH1
Species	Human
Gene ID (Human)	3269
Gene ID	3269
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-340 of human HRH1 (NP_000852.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P35367
Protein Weight	56kDa
Shipping	Ice bag