


Anti-Hrpt2 / Cdc73 Rabbit Polyclonal Antibody For Wb And Elisa - Q6P1J9

Item Number: CM0018395



Introduction

HRPT2 / CDC73 rabbit polyclonal antibody validated for Western blot (WB) and ELISA in human samples. Targets parafibromin, a tumor suppressor involved in transcription and Wnt signaling.

[Learn More](#)

Parameter	Value
Product Name	HRPT2 / CDC73 Rabbit pAb
Remarks/Alias	HYX; FIHP; HPTJT; HRPT1; HRPT2; C1orf28; HRPT2 / CDC73
Species	Human
GeneID (Human)	79577
GeneID	79577
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 435-505 of human HRPT2 / CDC73 (Q6P1J9)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6P1J9
Protein Weight	61kDa
Shipping	Ice bag