

Anti-Hrpt2/Cdc73 Rabbit Monoclonal Antibody For Wb And Elisa - Q6P1J9

Item Number: CM0004857

Introduction

Rabbit monoclonal antibody against HRPT2/CDC73 (Parafibromin). Validated for Western blot and ELISA in human, mouse, and rat. Tumor suppressor involved in transcription regulation and cell cycle progression.

[Learn More](#)

Parameter	Value
Product Name	HRPT2 / CDC73 Rabbit mAb
Remarks/Alias	HYX; FIHP; HPTJT; HRPT1; HRPT2; C1orf28; HRPT2 / CDC73
Species	Human
GeneID (Human)	79577
GeneID	79577
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-400 of human HRPT2 / CDC73 (NP_078805.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P1J9
Protein Weight	61 kDa
Shipping	Ice bag