

Anti-Hs1Bp3 Polyclonal Antibody For Wb, Elisa - Q53T59

Item Number: CM0025026

Introduction

Rabbit polyclonal antibody against human HS1BP3, suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Potential modulator of IL-2 signaling. Molecular weight 43 kDa.

[Learn More](#)

Parameter	Value
Product Name	HS1BP3 Rabbit pAb
Remarks/Alias	ETM2; HS1-BP3; HS1BP3
Species	Human
GeneID	64342
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human HS1BP3 (NP_071905.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q53T59
Protein Weight	43kDa
Shipping	Ice bag