

Anti-Hs3St3A1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y663

Item Number: CM0022546



Introduction

Rabbit polyclonal antibody against human HS3ST3A1, validated in WB and ELISA. Heparan sulfate 3-O-sulfotransferase involved in glycobiology and HSV-1 entry research.

[Learn More](#)

Parameter	Value
Product Name	HS3ST3A1 Rabbit pAb
Remarks/Alias	3OST3A1, 3-OST-3A, HS3ST3A1
Species	Human
Gene ID (Human)	9955
Gene ID	9955
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 44-140 of human HS3ST3A1 (NP_006033.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y663
Protein Weight	45kDa
Shipping	Ice bag