

Anti-Hsc70-Interacting Protein (Hip) Monoclonal Antibody For Wb - P50502

Item Number: CM0008487

Introduction

Anti-Hsc70-interacting protein (HIP) Rabbit Monoclonal Antibody for WB and ELISA. Recognizes human and mouse ST13/HIP, a cytoplasmic co-chaperone of Hsp70. 41 kDa protein. UniProt P50502.

[Learn More](#)

Parameter	Value
Product Name	Hsc70-interacting protein (HIP) Rabbit mAb
Remarks/Alias	HIP; HOP; P48; AAG2; SNC6; HSPABP; FAM10A1; FAM10A4; HSPABP1; PRO0786; Hsc70-interacting protein (HIP)
Species	Human
GeneID (Human)	6767
GeneID	6767
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 270-369 of human Hsc70-interacting protein (HIP) (P50502).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50502
Protein Weight	41kDa
Shipping	Ice bag