

Anti-Hsc70-Interacting Protein (Hip) Rabbit Pab - P50502

Item Number: CM0029145

Introduction

Rabbit polyclonal antibody recognizing human Hsc70-interacting protein (HIP/ST13) for Western blot (WB) and ELISA. Validated in human, mouse, rat. Suitable for chaperone research and protein folding studies.

[Learn More](#)

Parameter	Value
Product Name	Hsc70-interacting protein (HIP) Rabbit pAb
Remarks/Alias	HIP; HOP; P48; AAG2; SNC6; HSPABP; FAM10A1; FAM10A4; HSPABP1; PRO0786; Hsc70-interacting protein (HIP)
Species	Human
GeneID (human)	6767
GeneID	6767
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-369 of human Hsc70-interacting protein (HIP) (NP_003923.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50502
Protein Weight	41kDa
Shipping	Ice bag