

## Anti-Hsc70/Hspa8 Rabbit Monoclonal Antibody For Wb, Elisa - P11142

Item Number: CM0002936



### Introduction

High-quality rabbit monoclonal antibody against human Hsc70/HSPA8 (P11142). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for protein quality control, autophagy, and stress response research.

[Learn More](#)

| Parameter        | Value                                                                                                    |
|------------------|----------------------------------------------------------------------------------------------------------|
| Product Name     | Hsc70/HSPA8 Rabbit mAb                                                                                   |
| Remarks/Alias    | LAP1; HSC54; HSC70; HSC71; HSP71; HSP73; LAP-1; NIP71; HEL-33; HSPA10; HEL-S-72p; Hsc70/HSPA8            |
| Species          | Human                                                                                                    |
| GeneID (Human)   | 3312                                                                                                     |
| GeneID           | 3312                                                                                                     |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 547-646 of human Hsc70/HSPA8 (P11142) |
| Source           | Rabbit                                                                                                   |
| Category         | Monoclonal Antibody                                                                                      |
| Application      | WB, ELISA                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                        |
| SWISS            | P11142                                                                                                   |
| Protein Weight   | 71kDa                                                                                                    |
| Shipping         | Ice bag                                                                                                  |