

Anti-Hscb Rabbit Polyclonal Antibody For Wb And Elisa - Q8Iwl3

Item Number: CM0015659



Introduction

Rabbit polyclonal antibody against human HSCB (JAC1/HSC20), suitable for Western blotting and ELISA. Ideal for research on mitochondrial and cytoplasmic iron-sulfur cluster assembly, SDHB incorporation, and hematopoiesis.

[Learn More](#)

Parameter	Value
Product Name	HSCB Rabbit pAb
Remarks/Alias	JAC1; HSC20; SIDBA5; DNAJC20; HSCB
Species	Human
Gene ID (Human)	150274
Gene ID	150274
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-235 of human HSCB (NP_741999.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8IWL3
Protein Weight	27kDa
Shipping	Ice bag