

Anti-Hsd11B2 Rabbit Polyclonal Antibody For Wb, Elisa - P80365

Item Number: CM0027233

Introduction

Rabbit polyclonal antibody targeting human HSD11B2 (11-beta-hydroxysteroid dehydrogenase type 2). Validated for WB and ELISA, cross-reacts with rat. Suitable for glucocorticoid metabolism and mineralocorticoid receptor studies.

[Learn More](#)

Parameter	Value
Product Name	HSD11B2 Rabbit pAb
Remarks/Alias	AME; AME1; HSD2; HSD11K; SDR9C3; HSD11B2
Species	Human
Gene ID	3291 (Human)
Immunogen	Recombinant fusion protein containing amino acids 266-405 of human HSD11B2 (NP_000187.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P80365
Protein Weight	44kDa
Shipping	Ice bag