


Anti-Hsd17B1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P14061

Item Number: CM0018620



Introduction

HSD17B1 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA applications. Recognizes human, mouse, rat HSD17B1 (17 β -HSD1, 35kDa). Essential tool for studying estrogen biosynthesis and steroid metabolism.

[Learn More](#)

Parameter	Value
Product Name	HSD17B1 Rabbit pAb
Remarks/Alias	E2DH; HSD17; EDHB17; EDH17B2; SDR28C1; 17-beta-HSD; 20-alpha-HSD; HSD17B1
Species	Human
Gene ID (Human)	3292
Gene ID	3292
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-328 of human HSD17B1 (NP_000404.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14061
Protein Weight	35kDa
Shipping	Ice bag