

Anti-Hsd17B13 Rabbit Monoclonal Antibody For Wb, Elisa - Q7Z5P4

Item Number: CM0003110



Introduction

High-quality rabbit monoclonal antibody targeting human HSD17B13 (Gene ID: 345275). Validated for Western blot and ELISA. Ideal for hepatic lipid metabolism and liver disease research. SWISS: Q7Z5P4; Protein Weight: 34kDa.

[Learn More](#)

Parameter	Specification
Product Name	HSD17B13 Rabbit mAb
Remarks/Alias	FLDP; SCDR9; NIIL497; SDR16C3; HMFN0376; HSD17B13
Species	Human
GeneID (Human)	345275
GeneID	345275
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-211 of human HSD17B13 (NP_835236.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q7Z5P4
Protein Weight	34kDa
Shipping	Ice bag