

Anti-Hsd17B13 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q7Z5P4

Item Number: CM0027764



Introduction

Rabbit polyclonal antibody against HSD17B13, validated for WB, IF-P, IHC-P, ELISA. Detects human, mouse, rat HSD17B13; targets hepatic lipid metabolism research. Specific for lipid droplet and ER localization.

[Learn More](#)

Parameter	Value
Product Name	HSD17B13 Rabbit pAb
Remarks/Alias	FLDP; SCDR9; NIIL497; SDR16C3; HMFN0376; HSD17B13
Species	Human
GeneID (Human)	345275
GeneID	345275
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7Z5P4
Protein Weight	34kDa
Shipping	Ice bag