

Anti-Hsd17B13 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q7Z5P4

Item Number: CM0020362



Introduction

HSD17B13 Rabbit Polyclonal Antibody validated for WB, IF/ICC, IF-P, and ELISA. Detects human and rat HSD17B13, a ~34 kDa hepatic lipid metabolism enzyme. Ideal for liver research.

[Learn More](#)

Parameter	Value
Product Name	HSD17B13 Rabbit pAb
Remarks/Alias	FLDP, SCDR9, NIIL497, SDR16C3, HMFN0376, HSD17B13
Species	Human
GeneID (Human)	345275
GeneID	345275
Immunogen	Recombinant protein corresponding to a sequence within amino acids 34-300 of human HSD17B13 (NP_835236.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q7Z5P4
Protein Weight	34kDa
Shipping	Ice bag