

Anti-Hsd17B2 Monoclonal Antibody For Wb, Ihc-P, Elisa - P37059

Item Number: CM0012475

Introduction

Monoclonal antibody targeting human HSD17B2 (P37059), validated for WB, IHC-P, ELISA in human, mouse, rat. Ideal for studying steroid metabolism and endocrine regulation.

[Learn More](#)

Parameter	Value
Product Name	HSD17B2 Rabbit mAb
Remarks/Alias	HSD17; SDR9C2; EDH17B2
Species	Human
Gene ID (Human)	3294
Gene ID	3294
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human HSD17B2 (NP_002144.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37059
Protein Weight	43kDa
Shipping	Ice bag