

## Anti-Hsd17B3 Rabbit Polyclonal Antibody For Wb And Elisa - P37058

Item Number: CM0025034

### Introduction

High-quality rabbit polyclonal antibody to human HSD17B3, for WB and ELISA. Cross-reacts with mouse and rat. Useful for studying androgen biosynthesis and related diagnostic applications.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | HSD17B3 Rabbit pAb                                                                                                 |
| Remarks/Alias    | EDH17B3; SDR12C2; HSD17B3                                                                                          |
| Species          | Human                                                                                                              |
| Gene ID          | 3293                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human HSD17B3 (NP_000188.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | P37058                                                                                                             |
| Protein Weight   | 35kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |