

Anti-Hsd17B7 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P56937

Item Number: CM0022981

Introduction

Rabbit polyclonal antibody targeting human HSD17B7 (P56937), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse. Useful for steroid hormone metabolism and cholesterol biosynthesis research.

[Learn More](#)

Parameter	Value
Product Name	HSD17B7 Rabbit pAb
Remarks/Alias	PRAP; SDR37C1; HSD17B7
Species	Human
Gene ID (Human)	51478
Gene ID	51478
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-200 of human HSD17B7 (NP_057455.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P56937
Protein Weight	38kDa
Shipping	Ice bag