

Anti-Hsd17B8 Polyclonal Antibody For Wb, If/Icc, Elisa - Q92506

Item Number: CM0017222



Introduction

Rabbit polyclonal antibody targeting human HSD17B8 (Q92506), suited for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying mitochondrial fatty acid synthesis and steroid metabolism.

[Learn More](#)

Parameter	Value
Product Name	HSD17B8 Rabbit pAb
Remarks/Alias	KE6; FABG; HKE6; FABGL; RING2; H2-KE6; SDR30C1; D6S2245E; dJ1033B10.9; HSD17B8
Species	Human
GeneID (Human)	7923
GeneID	7923
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92506
Protein Weight	27kDa
Shipping	Ice bag