

Anti-Hsd3B1 Polyclonal Antibody For Wb, If-P, Elisa - P14060

Item Number: CM0026048

Introduction

Rabbit polyclonal antibody against human HSD3B1, validated for WB, IF-P, and ELISA. Detects endogenous HSD3B1 in human, mouse, and rat samples. Ideal for steroidogenesis and androgen metabolism research.

[Learn More](#)

Parameter	Value
Product Name	HSD3B1 Rabbit pAb
Remarks/Alias	HSD3B; HSDB3; HSDB3A; SDR11E1; 3BETAHSD; HSD3B1
Species	Human
GeneID (human)	3283
GeneID	3283
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human HSD3B1 (NP_000853.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14060
Protein Weight	42kDa
Shipping	Ice bag