

Anti-Hsd3B1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P14060

Item Number: CM0002211

Introduction

Rabbit monoclonal antibody against human HSD3B1, suitable for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Essential for steroidogenesis and hormone research.

[Learn More](#)

Parameter	Value
Product Name	HSD3B1 Rabbit mAb
Remarks/Alias	HSD3B; HSDB3; HSDB3A; SDR11E1; 3BETAHSD
Species	Human
GeneID (Human)	3283
GeneID	3283
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 250-350 of human HSD3B1 (P14060)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P14060
Protein Weight	42kDa
Shipping	Ice bag