

Anti-Hsd3B2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P26439

Item Number: CM0023541

Introduction

Rabbit polyclonal antibody against human HSD3B2, validated for Western blotting (WB), immunofluorescence (IF/ICC), and ELISA. Detects human and mouse HSD3B2. Ideal for research on steroidogenesis, hormone biosynthesis, and endocrine disorders.

[Learn More](#)

Parameter	Value
Product Name	HSD3B2 Rabbit pAb
Aliases	HSDB; HSD3B; SDR11E2; HSD3B2
Species	Human
Gene ID	3284
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-287 of human HSD3B2 (NP_000189.1)
Source	Rabbit
Category	Polyclonal Antibodies
Applications	WB, IF/ICC, ELISA
Cross-Reactivity	Human, Mouse
UniProt Accession	P26439
Molecular Weight	42kDa
Shipping	Ice bag