


Anti-Hsf1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q00613

Item Number: CM0004834



Introduction

Rabbit anti-HSF1 monoclonal antibody validated for WB, IHC-P, IP, ELISA. Reacts with human, mouse, rat. Ideal for heat shock response and transcription research. SWISS: Q00613.

[Learn More](#)

Parameter	Value
Product Name	HSF1 Rabbit mAb
Remarks/Alias	HSTF1; HSF1
Species	Human
GeneID	3297
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 429-530 of human HSF1 (NP_005517.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00613
Protein Weight	57kDa
Shipping	Ice bag