


Anti-Hsf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q00613

Item Number: CM0019642



Introduction

HSF1 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA in Human and Mouse. Stress-inducible transcription factor central to heat shock response. Suitable for cancer, stress biology, and HIV latency research.

[Learn More](#)

Parameter	Value
Product Name	HSF1 Rabbit pAb
Remarks/Alias	HSF1
Species	Human
GeneID (human)	3297
GeneID	3297
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 367-529 of human HSF1 (NP_005517.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q00613
Protein Weight	52kDa/57kDa
Shipping	Ice bag