


Anti-Hsf2 Rabbit Pab For Wb, Elisa - Q03933

Item Number: CM0021852



Introduction

High-quality rabbit polyclonal antibody targeting human HSF2 (Heat shock factor protein 2). Validated for WB and ELISA, reacts with human, mouse, and rat. Ideal for stress response and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	HSF2 Rabbit pAb
Remarks/Alias	HSF 2; HSTF 2; HSF2
Species	Human
GeneID (Human)	3298
GeneID	3298
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 39-215 of human HSF2 (Q03933).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03933
Protein Weight	60kDa
Shipping	Ice bag