


Anti-Hsl Rabbit Polyclonal Antibody For Wb, Elisa - Q05469

Item Number: CM0020663



Introduction

HSL rabbit polyclonal antibody for WB and ELISA applications, reactive with human, mouse, and rat HSL. Hormone-sensitive lipase (HSL) plays a key role in lipid metabolism by hydrolyzing triglycerides and cholesteryl esters. Ideal for metabolic studies.

[Learn More](#)

Parameter	Value
Product Name	HSL Rabbit pAb
Remarks/Alias	HSL; LHS; REH; AOMS4; FPLD6
Species	Human
Gene ID (Human)	3991
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 900-1000 of human LIPE (NP_005348.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05469
Protein Weight	117kDa
Shipping	Ice bag