


Anti-Hsl/Lipe Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q05469

Item Number: CM0007105



Introduction

Rabbit monoclonal antibody against human HSL/LIPE. Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse. Ideal for lipid metabolism, hormone-sensitive lipase, and steroidogenic research.

[Learn More](#)

Parameter	Value
Product Name	HSL/LIPE Rabbit mAb
Remarks/Alias	HSL; LHS; REH; AOMS4; FPLD6; HSL/LIPE
Species	Human
Gene ID (Human)	3991
Gene ID	3991
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 901-1000 of human HSL/LIPE (NP_005348.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q05469
Protein Weight	117kDa
Shipping	Ice bag