

Anti-Hsp105/Hsph1 Rabbit Monoclonal Antibody For Wb And Elisa - Q92598

Item Number: CM0009943



Introduction

Rabbit monoclonal antibody against human heat shock protein 105 kDa (HSP105/HSPH1), suitable for Western Blot and ELISA. Cross-reacts with mouse and rat. Ideal for stress response and protein folding research.

[Learn More](#)

Parameter	Value
Product Name	HSP105/HSPH1 Rabbit mAb
Remarks/Alias	HSP105; HSP105A; HSP105B; NY-CO-25; HSP105/HSPH1
Species	Human
GeneID (Human)	10808
GeneID	10808
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human HSP105/HSPH1 (Q92598)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92598
Protein Weight	97kDa
Shipping	Ice bag