

Anti-Hsp105/Hsph1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q92598

Item Number: CM0030439



Introduction

High-quality rabbit polyclonal antibody against human HSP105/HSPH1 (HSP105A/B, NY-CO-25). Validated for WB, IHC-P, ELISA with human, mouse, rat reactivity. Ideal for stress response and chaperone research.

[Learn More](#)

Parameter	Value
Product Name	HSP105/HSPH1 Rabbit pAb
Remarks/Alias	HSP105; HSP105A; HSP105B; NY-CO-25; HSP105/HSPH1
Species	Human
Gene ID (Human)	10808
Gene ID	10808
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92598
Protein Weight	97kDa
Shipping	Ice bag