


Anti-Hsp20/Hspb6 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - O14558

Item Number: CM0011233



Introduction

Rabbit monoclonal antibody targeting human HSP20/HSPB6 (O14558), validated for WB, IF-P, IHC-P, and ELISA. Suitable for research on molecular chaperone function, stress response, and cardioprotection.

[Learn More](#)

Parameter	Description
Product Name	HSP20/HSPB6 Rabbit mAb
Remarks/Alias	HEL55; Hsp20; PPP1R91; HSP20/HSPB6
Species	Human
GeneID (Human)	126393
GeneID	126393
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 70-150 of human HSP20/HSPB6 (O14558)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14558
Protein Weight	17kDa
Shipping	Ice bag